

DERMAQUAL

Designed to perform

INNOVATIVE DEVICES TO PROMOTE TREATMENT WITH GREAT RESULTS

FORZA XPERT MESO INJECTOR SYSTEM

The **FORZA XPERT** is a revolutionary aesthetic device that combines nano-penetration and electroporation technology as a single treatment. Dermaqual's Electroporation Microchip System has proven itself in professional, highly advanced therapies. The increased infusion of powerful vitamins into the mesoderm provides visible, natural results during aesthetic and dermatological procedures,





Micro-Penetration

The nano-needles cross the skin barrier to the mesoderm level, opening microchannels for the absorption of meso solutions, without causing skin irritation. highly advanced therapies. The increased infusion of powerful vitamins into the mesoderm provides visible, natural results during aesthetic and dermatological procedures.



Electroporation

Electroporation promotes the penetration of meso solutions, increases the absorption ratio, and helps deliver nutrients to subcutaneous skin tissue

ADVANCED RESULTS

The Multi-Channel Electroporation Microchip targets the epidermal layer, providing the best aesthetic results for clinical standards.











After two months of treatment, there was significantly reduced hair loss, and reactivated hair growth.

FORZA MMS 007

AUTOMATED MICRO NEEDLING & AUTO FUSION

The latest-technology, CE-mark-certified, automated, meso microneedling device employs a microchip with multiple pins to deliver nutrition and essentials to the skin, reaching the subcutaneous layer under high pressure from the pistol. The biggest advantage of the device is the auto-fusion system, which provides consistent application of the active ingredients in one shot. FORZA MMS 007 maintains continuous pressure in the syringe, allowing control button to economize the excess waste of meso-solutions and improving the results of classical mesotherapy treatments. The pressure from the bottom of the needle cartridge works like a very small dermo-jet, and in this way, multiple micro-injections are released without damaging the cutaneous. The microneedles are used to penetrate the skin, open the pores of the epidermal layer, and inject the skin boosters into the mesoderm.

Provides consatnt application of the active ingredients in one shot

FEATURES & ADVANTAGES

5 Speeds

Adjustable cartridge needle length: 0.25 - 2.5mm

Works wirelessly or with a power cord

Creation of thousands of microchannels in an effective, fast, and painless manner

Uninterrupted application of meso-active solutions

Combines active ingredients with safe, hygienic standards, and minimizes product waste

No cross-contamination (disposable needles)

Minimally invasive, no downtime

Suitable for all skin types

Short treatment time

DESIGNED FOR

Antiaging: Reduction of wrinkles and fine lines.

Rejuvenation: Improvement of skin elasticity and hydration levels.

Minimize the appearance of acne scars / pore refinement.

Improve skin tone, texture, and pigmentation.

Stretch marks and cellulite treatment.

Hair loss treatment/re-growth.

NEEDLE CARTRIDGES

12 NEEDLE

24 NEEDLE

......

The pressure from the bottom of the needle cartridge works like a very small dermo-jet, and in this way, multiple micro-injections are released without damaging the cutaneous.





FORZA MMS 600

AUTOMATED MICRO-NEEDLING

FORZA MMS 600 – One of the most advanced automated meso micro-needling and derma-stamp devices, offering high-speed oscillation, a powerful wireless motor, and an adjustable needle length. The internal motor causes the needles in the tip to oscillate at a rate of close to 110 revolutions per second. Utilizing Uniform Motion Technology, the derma pen FORZA MMS 600 penetrates the skin with a disposable needle cartridge, painlessly and efficiently creating thousands of micro-channels in a few seconds. When the needle tip makes contact with the skin, the needles create tiny punctures in the skin's surface, causing micro-trauma and inducing an immune response that subsequently creates new healthy collagen. The punctures also act as micro-channels that carry topical nutrients used during the treatment deep into skin lavers for an increased effect.

PACKAGE INCLUDES

Host / 1 pc Battery / 2 pcs Cord / 1 pc Charger / 1 pc Cover / 1 pc Tray / 1 pc

Adapter (110/220 V) general use / 1pc Needle cartridges (sterile, single use) /2 pcs – 36 pin

Product information guide / 1pc

Package box / 1pc

FEATURES & ADVANTAGES

Safer and more effective than traditional skin needling rollers for piercing the skin.

Various micro-needle modes, and uniform puncture-depth control across all treated areas. Administration of meso-active incredients.

Extremely high absorption of any active ingredients.

Microdermabrasion.

Multi-speed and adjustable needle depth (0.25 – 2.5 mm) significantly improves treatment. Utilizes sterile single-use surgical-grade disposable needle tips. No cross-contamination.

Shorter treatment time, with no side effects.

Painless procedure and greater comfort for the patient.

Visual and noticeable skin improvement in a few days from the first treatment.

Suitable for all skin types.

ONGOING & LONG-LASTING RESULTS

Superior effect on natural collagen and elastin reproduction.

Skin rejuvenation / anti-aging / fine lines and wrinkle reduction.

Reduces the appearance of acne scars, pimples, and blemishes.

Minimizes the appearance of large pores and blackheads.

Deep scar / pore refinement.

Enhance / nutrition skin treatment.

Removes stretch marks.

Hair loss treatment/re-growth.





DQ GOLD STAMP

NANO DELIVERY STAMP SYSTEM

DQ GOLD STAMP Nano Delivery Stamp System by DERMAQUAL is a new minimally invasive procedure involving superficial and controlled puncturing of the skin with fine miniature needles. DQ GOLD STAMP works by simply stamping the head full of safe and effective micro needles onto the skin. The micro-channel needles penetrate right into the dermis or the middle layer of the skin, this triggers a healing response that boosts production of natural collagen, the protein fibers that keep skin smooth and supple.

FEATURES & ADVANTAGES

20 Needles Titanium + 24K Gold Plated causes vertical penetration, creating infusion channels in the skin allowing solution to access deep into the epidermis maximizing absorption up to 98%.

Skin improvements is significantly enhanced by combining DQ GOLD STAMP with DERMAQUAL mesotherapy cocktalls according to the client's skin type or concern.

BENEFITS

Anti-Aging
Tightens skin
Improves skin texture and tone
Stimulates skin repair and cell generation
Smooths skin scarring
Reduces the appearance of pores
Reduces hair loss







SPECIFICATIONS

- Vial Material: Glass, sterile, color clear. Capacity: 6.5 ml / 0.22 fl.oz.
- Hypoallergenic Package Included: 1PC, sterile / BOX.
- . May be used on all skin types.

WHAT AREAS CAN BE TREATED WITH DQ GOLD STAMP?



Derma stamping is extremely versatile and can be used on all parts of the body including the face, neck, décolleté, scalp, arms, hands, legs, abdomen and back.

DQ GOLD STAMP is an effective way to combat hair loss.

The device creates micro-injury to the scalp, which is a catalyst for the body's natural healing response, that response instigates collagen production — which in turn stimulates hair growth. The process of derma stamping stimulates blood cells while improving new growth and the overall health of follicles.





DQ GOLD STAMD

DQ GOLD ROLLER

NANO DEL IVERY BOLLER SYSTEM

A new minimally invasive procedure involving superficial and controlled puncturing of the skin with fine miniature needles. The micro-channels that the treatment creates increases the rate of absorption and effectiveness of any subsequent products placed on the skin. A procedure that is also referred to as "collagen induction therapy," which aims to restore youth and smoothness to the skin by stimulating the natural production of collagen.

SPECIFICATIONS



- Vial Material: Glass, sterile, color clear. Capacity: 10 ml / 0.34 fl.oz.
- Hypoallergenic Package Included: 1PC, sterile / BOX.
- May be used on all skin types.



FEATURES

64 Needles Titanium + 24K Gold Plated allow solution to access deep into the epidermis through microscopic channels, maximizing absorption up to 98%.

DQ GOLD ROLLER 64 needles handpiece covers the largest treatment area for faster, more efficient treatments.

No ablation of the epidermis — no skin layers are removed, less chance of side effects, less pain and minimal healing time.

Improves product absorption. **DQ GOLD ROLLER** Nano Delivery Roller System can help mesotherapy penetrate skin more efficiently.



DQ GOLD ROLLER IS INTRODUCED FOR: Meso solution delivery, skin rejuvenation and tightening, pigmentation, acne scarring, stretch marks, and PRP use. DQ GOLD ROLLER Nano Delivery Roller System is the most gentle and effective skin treatment in mesotherapy.

WHAT AREAS CAN BE TREATED WITH DO GOLD BOLLER?

Microneedling is extremely versatile and can be used on all parts of the body including the face, neck, décolleté, arms, hands, legs, abdomen and back.



ICE ROLLER

PROFESSIONAL CRYO DEVICE

Is a professional Cryo Handheld Device with a water and gel filled rolling head. This roller uses an ice-cold barrel that rolls smoothly over the skin providing relaxation and relief instantly. DQ ICE ROLLER can be used to reduces pain, puffiness, redness and tighten skin. The cold therapy approaches facial rejuvenation and shortens time and truly making a "quicker lift."



INDICATIONS

It can be used before / after cosmetic procedures: injectables (mesotherapy, botox, filler).

laser, chemical peel, etc). It minimizes discomforts, irritation, and redness.

Speeds up recovery time of all cosmetic treatments.

 Calm your skin redness and inflammation - The detachable Cryo head holds the cold for up to 2 hours.

Effectively accelerate microcirculation on treatment area to clam your skin, reduce redness, fight puffy, fatique, dry eyes, dark circles, early fine lines.

Help tighten the jawline, minimize the appearance of large pores, fine lines and wrinkles.

Can be used on the body.

Relieves pain, tension, and discomfort associated with waxing or laser therapy.

Suitable for all skin types.

ONGOING & LONG-LASTING RESULTS

- 1. Detach the roller head, sanitize it with cotton, and 70% alcohol.
- 2. Store the roller head (in the zip bag) in the freezer for 4 hours or place it in the freezer for 15 minutes for a quick freeze.
- 3. Attach the freezed roller head.
- 4. Once cold, gently roll the device on desired areas such as face, neck, chest or body with a back and forth rolling motion approximately 2-5 minutes. Do not hold still on skin for more than 30 seconds. Apply to clean skin only.
- Wipe the roller head clean with 70% alcohol and place back in the freezer for next use.
- 6. Use as often as needed.







DERMAQUAL

MANUFACTURER AND MARKETING AUTHORIZATION HOLDER

DERMAQUAL 78/76 Ave. des Champs-Élysées 75008, Paris, FRANCE. Tel: 28 22 78 84 1(0) 33+ www.dermaqual.com info@dermaqual.com